


# Applications

## MecaTract




1) Fluorescence image of the beads near the surface of the hydrogel.  
2) Intensity profile of the fluorescent beads.

Michel Moussus' thesis (LTM)



HUVECS cells transfected with LifeAct Dye (Actin).  
Beads  $0.2 \mu\text{m}$ .


Courtesy of Alice Nicolas (LTM lab)



Intracellular stress pattern.  
HUVECS, 5 kPa.  
Soft Matter. Moussus et al, 2014


## MecaChips<sup>®</sup>

Brain Mecachips<sup>™</sup>




Mice neuronal cells growth on soft Laminin/ Poly-L-Lysine coated matrix.

Heart Mecachips<sup>™</sup>



iPSC-derived contractile cardio-myocyte on soft fibronectin coated matrix.


Lung Mecachips<sup>™</sup>



Human bronchial epithelial cells on soft collagen I coated matrix.


## $\mu$ pattern MecaChips<sup>®</sup>

Sub-cellular scale




REF52 cell line YFP paxillin. Dots  $3 \mu\text{m}$ .

Cellular scale



REF52 cell line  
 $\mu$ pattern  $80 \mu\text{m}$  by  $40 \mu\text{m}$ .


Tissue scale




a) stiff stripes  $10 \mu\text{m}$ / soft stripes  $90 \mu\text{m}$ .  
b) stiff stripes  $10 \mu\text{m}$ /soft stripes  $40 \mu\text{m}$ .

Cell&Soft culture plates  
are dedicated to

Cell culture  
& Cell-based assays



Drug discovery



Stem cells



Microscopy &  
Time-lapse

## Contact

Cell&Soft SAS

c/o CEA/CNRS/LTM

17 Avenue des Martyrs

Bât 44

38054 Grenoble Cedex 9

[contact@cellandsoft.com](mailto:contact@cellandsoft.com)

+33 (0)4 38 78 48 10

[www.cellandsoft.com](http://www.cellandsoft.com)

Cell&Soft SAS property 2021

# Cell&Soft

*Let's reinvent cell culture*

## Product Catalog

*Soft culture plates  
for pleased cells*

Robust  
Repeatability

Reliable  
Reproducibility

Ready  
to use  
Time saving

Easy  
to use  
Use it as a  
standard plate



Soft plates for cell culture

# MecaTract

Cell-contractile forces generated by the actomyosin cytoskeleton and transmitted to the extracellular matrix (ECM) drive cell adhesion, spreading and migration. These forces are known to be critical during embryo morphogenesis, wound healing, immune response as well as pathological processes such as cancer metastasis. Traction force microscopy (TFM) is a recognized experimental technique that measures the surface forces, also termed as tractions, that cells exert on a given substrate. It relies on the computational analysis of the direction and magnitude of elastic substrate deformations to reconstruct cell-generated traction forces. These deformations can be tracked and quantified by recording the displacement of fluorescent beads already embedded in the substrate, as a result of the mechanical stress induced by an adherent cell.

## Applications

### Traction Force Microscopy

#### Fluorescent beads

Size: 0.2  $\mu\text{m}$   
Fluorophore: Dark red  
Wavelength (Exit. / Emiss.): 660/680 nm  
Cell seeding surface: 6.15  $\text{cm}^2$   
(Specific beads on request)

ONE FORMAT: PD35

#### BENEFITS

- Robust and reliable
- Ready and easy to use
- No biological risk (synthetic matrix)
- Glass bottom: ideal for microscopy
- Stiffness and coating are decoupled

**Polyacrylamide gel with fluorescent carboxylate-modified microspheres uniformly dispersed inside**

Soft plates for cell culture


# MecaChips<sup>®</sup>

*In vivo*, cells lay in soft tissues with distinct physical properties. Rigidity plays a major role in a myriad of cellular mechanisms, such as carcinogenesis and metastasis formation, as well as stem cells differentiation and drug effectiveness. Mecachips<sup>®</sup> soft and flat matrices are new and physiological solutions for *in vitro* cell culture. They mimic the soft mechanical features of all human or animal tissues, thus preserving the cells *in vivo* characteristics.

## Applications

**Areas:** basic research, stem cells, oncology, neurology, cardiology, compatible with HCS/HTS platforms.

**Tools:** cell biology, molecular biology, biochemistry.




#### Hydrogel + ECM

#### BENEFITS

- Scalable to screening workflow (no topography)
- No biological risk (synthetic matrix)
- Glass bottom: ideal for microscopy
- Stiffness and coating are decoupled
- Compatible with standard analysis

#### CHOOSE YOUR FORMAT\*



\*(P384 on request)

#### CHOOSE YOUR ORGAN OF INTEREST



## Mecachips<sup>®</sup> Starter Kit

Get a pack of 3 plates to familiarize yourself with our substrates:

- 1 Select 3 organs of interest
- 2 Choose your coating for each of them
- 3 Choose between 3 formats: PD35, P6, P96

STORAGE  
Temperature: +4°C  
Shelf life: 3 months

#### CHOOSE YOUR COATING

- Culture-dedicated surface chemistry:
- Vitronectin (human, recombinant truncated)
  - Fibronectin (human, plasma)
  - Collagen I (rat, tail)
  - Laminin (mouse, EHS sarcoma)
  - Poly-Ornithine/ Laminin
  - Poly-Ornithine

**To be representative, *in vitro* cell devices should selectively provide culture conditions as close as possible to the mechanical microenvironment of targeted tissues.**

Soft plates for cell culture

# $\mu$ pattern MecaChips<sup>®</sup>

*In vivo* tissues are soft, elastic and mechanically textured. Living tissues rigidity properties present microscale variations that can play a crucial role in cell response (in a same tissue, rigidity varies from Pa/mm to kPa/ $\mu\text{m}$ ). Relying on unique patented technologies and know-how derived from the microelectronic field, the mechanical properties of Mecachips<sup>®</sup> matrices can be finely tuned up to the  $\mu\text{m}$  scale to replicate such variations.

## Applications

**Areas:** basic research, stem cells, oncology, neurology, cardiology.

**Tools:** adhesion, cell shape standardization, cell confinement, cell migration.

